

Features

- UL® 1069 Listed
- Oxygen safe
- Fully sealed for easy cleaning
- Flame-retardant (UL® 94V-0) high-impact plastic
- Long-lasting switches surpass UL® 1069 requirements
- 8-foot cord with DIN plug and security clip
- Provides pillow speaker intercom and light control (when used with appropriate modules)
- Special order: analog volume control, 10-foot cord, earphone jack

Description

The SF41ZKL Pillow Speaker is used with Tek-CARE400 IR411-series and IR412-series patient stations, Tek-CARE400 P5 IR411P5 and IR412P5 patient stations, and Tek-CARE300II/III IR319-series and IR320-series patient stations.

The pillow speaker provides call placement to the master station, as well as remote control for hospital-grade TVs. The pillow speaker is programmed by default for Zenith Digital TVs, however, the pillow speaker can be reprogrammed to match the code set that is required by other hospital-grade TVs.

When used in conjunction with the appropriate TekTone pillow speaker modules, the pillow speaker also supplies two-way, hands-free communication with the master station, and automatic silencing of TV audio while the intercom is in use.

Push buttons are provided to call the nurse, plus TV power, channel selection up/down, Closed Captioning, mute and two light controls. Buttons are also provided for volume up/down (analog volume controls are available by special order) and for direct channel access (in some models).

Listings



UL® 1069 Listed and cUL® Listed to CSA C22.2 No. 205; UL® File Number E71599.



SF41ZKL Pillow Speaker

Required Components

Hospital-grade TV

for Tek-CARE®300II/III Nurse Call Systems:
IR319- or IR320-series Patient Station
PM311C or PM312C Pillow Speaker Intercom Module
PM321B or PM322B Light Control Module

for Tek-CARE®400 Nurse Call Systems:
IR411- or IR412-series Patient Station
PM421 or PM422 Pillow Speaker & Bed Interface Module

for Tek-CARE®400 P5 Nurse Call Systems:
IR411P5 or IR412P5 Patient Station
PM421P5 or PM422P5 Pillow Speaker Module