

### Features

- UL® 1069 Listed
- Oxygen safe
- Fully sealed for easy cleaning
- Flame-retardant (UL® 94V-0) high-impact plastic
- Long-lasting switches surpass UL® 1069 requirements
- 8-foot cord with 8-PIN 8P8C plug and security clip
- Provides pillow speaker intercom and light control (when used with appropriate modules)
- Special order: analog volume control, 10-foot cord, earphone jack

### Description

The SF31ZDL Pillow Speaker is used with Tek-CARE300III IR311-series and IR312-series patient stations. The pillow speaker provides call placement to the master station, as well as remote control for hospital-grade digital TVs (Zenith is standard; RCA and Philips are available by special order). When used in conjunction with the appropriate TekTone® pillow speaker modules, the pillow speaker also supplies two-way, hands-free communication with the master station, and automatic silencing of TV audio while the intercom is in use.

Push buttons are provided to call the nurse, plus TV power, channel selection up/down, closed captioning, mute and two light controls. Buttons are also provided for volume up/down (analog volume controls are available by special order), and for direct channel access (in some models).

### Listings



UL® 1069 listed; UL® File Number E71599.



SF31ZDL Pillow Speaker

### Required Components

Hospital-grade Zenith TV  
IR311- or IR312-series Patient Station  
PM311C or PM312C Pillow Speaker Intercom Module  
PM321B or PM322B Light Control Module