

# TekTone® Tek-CARE® Appliance Server Integration Layout

(See installation manuals for complete wiring information.)

**Tek-CARE500**  
Wireless Nurse Call

**Tek-CARE570**  
Wireless Nurse Call

**Tek-CARE400 GEN3**  
Nurse Call

**Tek-CARE120/160**  
Nurse Call

**Tek-CARE750**  
Fall Detection and Monitoring

**Tek-CARE NC110, NC150, NC200**  
Nurse Call

**Tek-CARE700**  
Wander Management

The Tek-CARE Appliance Server can simultaneously connect to any TekTone nurse call system, plus a variety of foreign systems—such as fire alarm and security systems. The Tek-CARE Appliance Server displays the calls, events and messages from all connected systems at once.

**NC475 Tek-CARE Appliance Server**

Any computer networked with the Tek-CARE Appliance Server can create customizable reports on all activity. Plus, networked PCs can view all system events and messages as they occur via TekTone's Tek-CARE Remote Event Monitor.

**Tek-ALERT**  
Integration Manager

foreign systems

Facility LAN

**Tek-PAGING**  
Radio Pocket Paging

**Tek-CARE Remote Event Monitor**

**Reporting System and Config Tool Live**

**Media Gateway License**

**Email Output Software**

**Staff App**  
iPhone, iPad, Android

**Tek-CARE TV Monitors**